


# BECO Validation Guide

**BECO**PAD P



# BECO Depth Filter Media

## **BECOPAD P**

UH, SBr  
11/2010  
Page 2 of 18

---

### Table of Contents

<b>1. Introduction</b>	page 3
<b>2. General Characteristics of BECOPAD P Range Depth Filter Media</b>	page 4
2.1 Overview on <b>BECOPAD P</b> types	
2.2 Filtration mechanisms	page 5
2.3 Physical Data	page 6
<b>3. Extractable Tests</b>	page 7
3.1. Conductivity	page 8
3.2. Extractable Cations (Ca, Mg, Fe, Al, Cr, Cu, Zn, Cd, Mn, As, Pb) in deionized water and 40 % ethanol	page 9
3.3. Total Extractables according to the guideline of FDA CFR 21 Chapter 177.2260	page 11
3.4. Investigation according to LFGB and recommendation XXXVI/1 of the BfR	
3.5. Total organic Carbon	page 12
3.6. Oxidizable Substances	page 13
<b>4. Endotoxins</b>	page 14
<b>5. Bacterial Challenge Testing</b>	page 15
<b>6. Chemical Data</b>	page 16
<b>7. Biological Safety</b>	page 17
<b>8. Appendix</b>	page 18
8.1 Appendix I (14 pages)	
8.2 Appendix II (7 pages)	

# BECO Depth Filter Media

## **BECOPAD P**

UH, SBr  
11/2010  
Page 3 of 18

---

### 1. Introduction

---

Depth filtration is regularly used for the separation of cells, cell fragments, colloids, precipitates, and other fermentation process components. Also, this technology is trend-setting within the filtration of media containing particles or for removal of bioburden. A further cost-effective application is the protection of final sterile filters from premature blocking using sterilizing depth filter sheets.

Clarifying depth filtration is a very important process step in the production of biological products, as it has a direct influence on the yield, product characteristics and reproducibility of the final product. Depth filtration therefore contributes to the economic efficiency and safety of the overall process.

We are committed to providing consistently high product quality in order to meet the quality requirements of our customers. As part of the manufacturing process, the depth filter medium helps to ensure safe, reproducible, and economical results for our customers.

We meet this responsibility by ensuring that our products comply with national and international quality standards such as the requirements of recommendation XXXVI/1 regarding the LFGB (German Food, Commodity and Feed Act) by the BfR (Federal Institute of Risk Assessment), and the test criteria of FDA directive CFR 21 § 177.2260 and the hygiene guidelines according to HACCP. As early as 1993, BEGEROW was certified according to DIN EN ISO 9001 and has been complying with DIN EN ISO 9001:2000 since 2002. Moreover, we make ourselves available for regular audits through external, independent institutes, and customers.

This **BECOPAD P** Validation Guide summarizes the extensive testing and evaluation program prepared by BEGEROW to meet the critical demands of the Pharmaceutical Industry.

Validation is defined as “Establishing documented evidence which provides a high degree of assurance that a specific process will consistently produce a product meeting its predetermined specifications and quality Attributes“ (Guideline For Industry: Sterile Drug Products Produced by Aseptic Processing - Current Good Manufacturing Practice (US Department of Health and Human Services, Food and Drug Administration, September 2004).

This report contains validation data applicable to BEGEROW **BECOPAD P** range depth filter media which is used in different formats and configurations such as:

- BECO Flat Sheets for use in filter presses
- BECODISC Stacked Disc Cartridges used in stacked disc cartridge housings
- Disposable BECO MiniCaps and BECO Capsules

This validation guide has been compiled for the users of BEGEROW **BECOPAD P** range depth filter media as a basis and support for their own validation procedures.

# BECO Depth Filter Media

## **BECOPAD P**

UH, SBr  
11/2010  
Page 4 of 18

---

## 2. General Characteristics of **BECOPAD P** Range Depth Filter Media

---

### 2.1 Overview of **BECOPAD P** types

The high-purity **BECOPAD P** range is a mineral-free depth filter medium produced to date with special, high-purity cellulose using the bepure process.

This validation report describes the following depth filter media grades, listed with increasing permeability:

<b>BECOPAD P</b>
<b>BECOPAD P 120</b>
<b>BECOPAD P 170</b>
<b>BECOPAD P 220</b>
<b>BECOPAD P 270</b>
<b>BECOPAD P 350</b>
<b>BECOPAD P 450</b>
<b>BECOPAD P 550</b>

All of these depth filter media grades are manufactured under special production conditions that guarantee the highest purity possible.

Through the use of special production methods, **BECOPAD P** depth filter media can be distinguished by very low release of extractables. Furthermore **BECOPAD P** ranging from **BECOPAD P 120** through **BECOPAD P 550** have been optimized for low endotoxin levels.

# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 5 of 18

---

## 2. General Characteristics of BECOPAD P Range Depth Filter Media

---

### 2.2 Filtration mechanisms

#### General filtration mechanisms:

The depth filter medium offers mechanical and adsorptive mechanisms. Particles are mechanically retained in the asymmetric cavity structure of the depth filter medium (narrowing pore structure towards the filtrate side). Particles that are significantly smaller than the pores of the depth filter medium are bound through adsorption. This is due to the positive zeta potential (electrokinetic potential), which retains negatively charged particles through adsorption.

#### Filtration mechanisms of BECOPAD P:

In addition to mechanical separation of particles by their size, **BECOPAD P** offers additional specific separation based on its cationic charge. The key factor is the positive zeta potential (electrokinetic potential), which retains negatively charged particles through adsorption.

When product flows through the cationic **BECOPAD P**, negatively charged particles accumulate at the positively charged internal surfaces. Several successive boundary layers with different charge densities are formed, the effect of which becomes weaker with increasing distance from the original charge.

The cationic **BECOPAD P** is recommended for specific reduction of negatively charged particles, for example endotoxins.

#### Ingredients or material of construction:

All components of the **BECOPAD P** depth filter media are FDA listed for food contact use in the Code of Federal Regulations (CFR), Title 21, 177.2260.

#### Components:

- Highly pure Cellulose
- Wet strength agents

# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 6 of 18

### 2. General Characteristics of BECOPAD P Range Depth Filter Media


#### 2.3 Physical data

##### Overview physical data

Type	Article No.	Nominal filtration rate [µm]	Thickness [mm]	Ash content [%]	Bursting strength, wet [kPa]	Water throughput at $\Delta p = 100 \text{ kPa}$ [ $\text{l m}^{-2} \text{ min}^{-1}$ ]	Endotoxin content <sup>(1)</sup> [EU/ml]
<b>BECOPAD P 120</b>	Q1112	0.1 – 0.3	3.9	< 1	> 150	55	< 0.025
<b>BECOPAD P 170</b>	Q1117	0.2 – 0.4	3.9	< 1	> 150	80	< 0.025
<b>BECOPAD P 220</b>	Q1122	0.3 – 0.5	3.9	< 1	> 150	100	< 0.025
<b>BECOPAD P 270</b>	Q1127	0.5 – 0.7	3.9	< 1	> 150	135	< 0.025
<b>BECOPAD P 350</b>	Q1135	0.7 – 1.0	3.9	< 1	> 150	160	< 0.025
<b>BECOPAD P 450</b>	Q1145	1.0 – 2.0	3.9	< 1	> 150	300	< 0.025
<b>BECOPAD P 550</b>	Q1155	2.0 – 3.0	3.9	< 1	> 150	700	< 0.025

The endotoxin content is determined after rinsing with 25 l/m<sup>2</sup> of WFI.

#### Water Flow Rates for BECOPAD P



## BECO Depth Filter Media **BECOPAD P**

UH, SBr  
11/2010  
Page 7 of 18

---

### 3. Extractable Tests

---

#### **Method**

In the filtration of biotech and pharmaceutical products, it is essential that product composition is not changed by filtration. An appropriate rinsing procedure after sterilization was therefore used to remove any extractable substances. The extraction media refer to the two most characteristic process conditions deionized water and additionally for the ion extractions 40 % ethanol.

For the extraction curves discussed in the following chapter, a filter with an effective filter area of 131 cm<sup>2</sup> as used. The flow rate was adjusted to 500 l/m<sup>2</sup>/h.

# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 8 of 18

### 3. Extractable Tests

#### 3.1 Conductivity in deionized water

##### Method


Conductivity is measured by a calibrated conductivity test device

##### Results

##### Conductivity Shift (in deionized water)

	Rinse volume [l/m <sup>2</sup> ]				
	After 3.8	After 15.3	After 22.9	After 53.4	After 106.9
<b>BECOPAD P 120</b>	67	7	4	2	1
<b>BECOPAD P 170</b>	57	5	3	2	2
<b>BECOPAD P 220</b>	65	7	4	2	2
<b>BECOPAD P 270</b>	93	15	9	4	2
<b>BECOPAD P 350</b>	49	7	3	2	1
<b>BECOPAD P 450</b>	84	32	10	4	2
<b>BECOPAD P 550</b>	104	29	7	3	2

#### Conductivity in $\mu\text{S}/\text{cm}$ at different rinsing volumes



The conductivity of the extract was  $< 20 \mu\text{S}/\text{cm}$  after a rinsing volume of  $22.9 \text{ l}/\text{m}^2$  and can be lowered further by additional flushing.



# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 9 of 18

### 3. Extractable Tests

#### 3.2 Extractable Cations (Ca, Mg, Fe, Al, Cr, Cu, Zn, Cd, Mn, As, Pb) in deionized water and 40 % Ethanol

##### Method

For the extraction curves discussed in the following chapter, a filter with an effective filter area of 131cm<sup>2</sup> was used. The flow rate was adjusted to 500 l/m<sup>2</sup> h.

Samples were taken after a extraction volume of:

- 0 – 1.9 l/m<sup>2</sup>
- After 10 l/m<sup>2</sup>
- After 25 l/m<sup>2</sup>
- After 50 l/m<sup>2</sup>

All **BECOPAD P** depth filters are have been rinsed according to the recommendation with 25l/m<sup>2</sup> with deionized water.

All cations are measured by AAS (atomic adsorption spectroscopy) by flame or graphite tube technique.

##### Results

##### Typical ion contents for **BECOPAD P**

##### 0 – 1.9 l/m<sup>2</sup>

	Deionized water*	Ethanol**
Calcium	< 1 ppm	< 1 ppm
Magnesium	< 0.5 ppm	< 0,5 ppm
Iron	< 5 ppb	< 5 ppb
Aluminum	< 5 ppb	< 5 ppb
Chromium	< 5 ppb	< 5 ppb
Copper	< 5 ppb	< 5 ppb
Zinc	< 5 ppb	< 5 ppb
Cadmium	< 5 ppb	< 5 ppb
Manganese	< 5 ppb	< 5 ppb
Arsenic	< 5 ppb	< 5 ppb
Lead	< 5 ppb	< 5 ppb

##### After 10 l/m<sup>2</sup>

	Deionized water*	Ethanol**
Calcium	< 25 ppb	< 25 ppb
Magnesium	< 25 ppb	< 25 ppb
Iron	< 5 ppb	< 5 ppb
Aluminum	< 5 ppb	< 5 ppb
Chromium	< 5 ppb	< 5 ppb
Copper	< 5 ppb	< 5 ppb
Zinc	< 5 ppb	< 5 ppb
Cadmium	< 5 ppb	< 5 ppb
Manganese	< 5 ppb	< 5 ppb
Arsenic	< 5 ppb	< 5 ppb
Lead	< 5 ppb	< 5 ppb

## BECO Depth Filter Media

### BECOPAD P

#### After 25 l/m<sup>2</sup>

	Deionized water*	Ethanol**
Calcium	< 25 ppb	< 25 ppb
Magnesium	< 25 ppb	< 25 ppb
Iron	< 5 ppb	< 5 ppb
Aluminum	< 5 ppb	< 5 ppb
Chromium	< 5 ppb	< 5 ppb
Copper	< 5 ppb	< 5 ppb
Zinc	< 5 ppb	< 5 ppb
Cadmium	< 5 ppb	< 5 ppb
Manganese	< 5 ppb	< 5 ppb
Arsenic	< 5 ppb	< 5 ppb
Lead	< 5 ppb	< 5 ppb

#### After 50 l/m<sup>2</sup>

	Deionized water*	Ethanol**
Calcium	< 25 ppb	< 25 ppb
Magnesium	< 25 ppb	< 25 ppb
Iron	< 5 ppb	< 5 ppb
Aluminum	< 5 ppb	< 5 ppb
Chromium	< 5 ppb	< 5 ppb
Copper	< 5 ppb	< 5 ppb
Zinc	< 5 ppb	< 5 ppb
Cadmium	< 5 ppb	< 5 ppb
Manganese	< 5 ppb	< 5 ppb
Arsenic	< 5 ppb	< 5 ppb
Lead	< 5 ppb	< 5 ppb

#### Conclusion

All evaluated ions are extracted in on a very low level. This apply for deionized water as well as for 40 % ethanol.

# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 11 of 18

### 3. Extractable Tests

#### 3.3 Total Extractables according to the guideline of FDA CFR 21 Chapter 177.2260

##### Method

Under Title 21 «Food and Drugs» (2), Chapter 177 the Code of Federal Regulations of the FDA (Food and Drug Administration) is concerned with indirect food additives: polymers. Subpart C. Substances for use only as components of articles intended for repeated use.

The 177.2260 refers to resin-bonded filters and states the limits for extractables in different extraction media and under different extraction conditions. The following limitations are stated here.

Extraction Solvent	Extraction Conditions	CFR Limits
deionized Water	2 h/boiling heat	< 4 % dry residue
50 % ethanol	2 h/60 °C	< 4 % dry residue
8 % ethanol	2 h/80 °C	< 4 % dry residue
5 % acetic acid	2 h/90 °C	< 4 % dry residue
n-hexane	2 h/boiling heat	< 0.5 % dry residue

Although these regulations are not primarily intended for pharmaceutical products, they provide additional supporting data for the suitability of filters for pharmaceutical applications.

##### Results and Conclusion

The level of extractables of all BEGEROW depth filter media was significantly below the CFR limits and therefore meets the requirements of 21 CFR 177.2260.

For a detailed report please refer to Appendix I.

#### 3.4 Investigation according to LFGB and recommendation XXXVI/1 of the BfR

##### Method

**BECOPAD P** meet the requirements of LFGB (German Food, Commodity, and Feed Act), and the Recommendation XXXVI/1 issued by BfR (German Federal Institute of Risk Assessment).

This test report includes:

Hot water extracts for the determination of:

Formaldehyde	< 0.3 mg/g
Heavy metals:	
Lead	< 0.3 µg/g
Cadmium	< 0.5 µg/g
Mercury	< 0.3 µg/g
Chromium	not detectable
Total Dry Residue	< 5 mg/g
Total Nitrogen	< 3 mg/g dry substance

Sensory tests

Germ Growth Inhibition Tests

Detection of the migration of Dichloropropanol	< 2 µg/g
Monochloropropandiol	< 12 µg/g

##### Results and Conclusion

**BECOPAD P** depth filters pass all tests according the requirements LFGB.

A detailed report and further validation information are available by contacting BEGEROW.

# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
 11/2010  
 Page 12 of 18

### 3. Extractables Test

#### 3.5 Total Organic Carbon


##### Method

To determine the quantity of organic extractables in BECOPAD the depth filter media were tested for total organic carbon (TOC).

Four BECOPAD P depth filter types have been tested exemplary for the whole range. The filters with an area of 131 cm<sup>2</sup> m<sup>2</sup> have been flushed with deionized water at a flow rate of 500 L / m<sup>2</sup> h. Samples of flushing water were taken after 1 L, 25 L, 50 L and 100 L and analyzed for TOC contents.

##### Result

TOC values in flushing water in ppb





# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 14 of 18

### 4. Endotoxins

In a medical sense, pyrogen means fever pathogen. Examples are Endotoxins of gram-negative bacteria. These are cell wall constituents that consist of a lipid and a polysaccharide component and are embedded in the outer membrane of gram-negative microbes. Quantitative proof of Endotoxins can be established using the LAL test (Limulus Amoebocyte Lysate) which is also used for the routine examination of the **BECOPAD P** depth filters.

#### Methods

The Endotoxin content of the samples is specified in EU/ml (Endotoxin unit) testing with the kinetic turbidimetric LAL method. The sample is incubated with the LAL reagent and the rate of the increase in turbidity during the time is taken. It is measured spectrophotometrically and compared to a standard curve. Prior to the measurement of the Endotoxin content, depth filters are rinsed with 25 l/m<sup>2</sup> WFI. The tests are performed by a accredited laboratory according to the guidelines of Ph.Eur. 2.6.14 [harmonized testing].

#### Results

Type	Article No.	Endotoxin content <sup>(1)</sup> [EU/ml]
<b>BECOPAD P 120</b>	Q1112	< 0.025
<b>BECOPAD P 170</b>	Q1117	< 0.025
<b>BECOPAD P 220</b>	Q1122	< 0.025
<b>BECOPAD P 270</b>	Q1127	< 0.025
<b>BECOPAD P 350</b>	Q1135	< 0.025
<b>BECOPAD P 450</b>	Q1145	< 0.025
<b>BECOPAD P 550</b>	Q1155	< 0.025

<sup>(1)</sup> The endotoxin content is determined after rinsing with 25 l/m<sup>2</sup> of WFI.

#### Conclusion

Endotoxin limit values at < 0.025 EU/ml are defined for the **BECOPAD P** range.

# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 15 of 18

### 5. Bacterial Challenge Testing

Fine depth filters have the potential to reduce the bioburden load in a solution. To evaluate their bioburden reduction the depth filter is challenged with a defined bacteria suspension and the titer reduction in the filtrate is determined.

The reduction capability of depth filters is influenced by several factors such as the physical and chemical properties of the fluid, the process parameters and the type of challenge organism.

#### Result

Log. Reduction values for BECOPAD P depth filter media.

Type	Nominal retention rate [ $\mu\text{m}$ ]	Log Reduction Value [LRV]
BECOPAD P 120	0.1 – 0.3	8.0 <sup>1</sup>
BECOPAD P 170	0.2 – 0.4	6.7 <sup>1</sup>
BECOPAD P 220	0.3 – 0.5	5.9 <sup>2</sup>
BECOPAD P 270	0.5 – 0.7	5.1 <sup>2</sup>
BECOPAD P 350	0.7 – 1.0	3.5 <sup>2</sup>

<sup>1</sup> Challenged with *Brevundimonas diminuta* (ATCC 19146)

<sup>2</sup> Challenged with *Serratia marcescens* (ATCC 14765)

# BECO Depth Filter Media

## BECOPAD P

UH, SBr  
11/2010  
Page 16 of 18

### 6. Chemical Data

The **BECOPAD P** depth filters are tested by putting them under the written temperature for at maximum 168 h in the chemical compound. If the time is below 168 h, this means that the depth filter was monitored as resistant at the maximum given value of time.

#### Chemical resistance of **BECOPAD P** depth filter media

Chemical compound	Max. tested temperature/ Contact time	Mechanical resistance	Chemical compound	Max. tested temperature/ Contact time	Mechanical resistance
<b>Caustic:</b>			<b>Organic solvents:</b>		
ammonia solution	25% 20 °C/168 h	X	acetone	20 °C/168 h	X
potassium hydroxide	30% 20 °C/ 48 h	(X)	butanol	20 °C/168 h	X
sodium hydroxide	30% 20 °C/ 24 h	-	cyclohexane	20 °C/168 h	X
	5% 40 °C/ 4 h	X	dimethyl sulphide	20 °C/168 h	X
	2% 40 °C/ 4 h	X	ethanol	20 °C/168 h	X
	1% 40 °C/ 4 h	X	ethylene glycol	20 °C/168 h	X
	0.5% 40 °C/ 4 h	X	ethyl methyl ketone	20 °C/168 h	X
			isopropanol	20 °C/168 h	X
<b>Acids:</b>			methanol	20 °C/168 h	X
acetic acid	25% 20 °C/168 h	X	N,N dimethyl formamide	20 °C/168 h	X
peracetic acid	0.1% 20 °C/168 h	X	n-hexane	20 °C/168 h	X
	0.2% 20 °C/ 168 h	X	tetrachloroethylene	20 °C/168 h	X
	0.5% 20 °C/ 168 h	X	toluene	20 °C/168 h	X
nitric acid	25% 20 °C/ 48 h	X	triethanolamine	20 °C/168 h	X
hydrochloric acid	25% 20 °C/168 h	X	xylene	20 °C/168 h	X
sulphuric acid	25% 20 °C/ 48 h	X			
citric acid	25% 20 °C/168 h	X	<b>Aqueous solutions:</b>		
			iron trichloride	25% 20 °C/168 h	X
			sodium hypochlorite free chlorine	12% 20 °C/168 h	X
			hydrogen peroxide	10% 20 °C/ 72 h	X

X = resistant

(X) = limited resistance

- = not resistant



## BECO Depth Filter Media **BECOPAD P**

UH, SBr  
11/2010  
Page 17 of 18

---

### 7. Biological Safety

---

**BECOPAD P** depth filter media have been evaluated for their biological properties. According to the following guidelines the depth filter media have been examined for cytotoxicity, systemic cytotoxicity, irritation and tissue compatibility:

- Cytotoxicity in accordance to DIN EN ISO 10993-5:1999
- USP Plastic Class VI test in accordance to USP 32 <88>

#### **Method**

In the Cytotoxicity Test extracts of **BECOPAD P** depth filter media are investigated for any growth inhibition of L929 mouse fibroblasts. All extracts of **BECOPAD P** depth filter media did not show any growth inhibition effect.

Further the **BECOPAD P** depth filter media have been tested according to USP Plastic Class VI test in systemic injection test, intracutaneous reactivity test and implantation test for their biological safety in test animals. In all test animals no effects were observed.

#### **Result**

Based on the test results described above it can be stated that **BECOPAD P** depth filter media have no cytotoxic potential and comply with the requirements of USP Plastic Class VI test in accordance to USP 32 <88>.

For details see Appendix II.

# **BECO Depth Filter Media**

## **BECOPAD P**

UH, SBr  
11/2010  
Page 18 of 18

---

### **8. Appendix**

---

#### **8.1 Appendix I (14 pages)**

Test reports "Total extractables acc. to the guideline of FDA CFR 21 Chapter 177.2260"  
Test reports acc. to LFBG

#### **8.2 Appendix II (7 pages)**

Test report "Cytotoxicity in accordance to DIN EN ISO 10993-5:1999"  
Test report "USP Plastic Class test in accordance to USP 32 <88>"