

## Depth Filtration

### BECODISC® B30C

#### Activated Carbon Stacked Disc Cartridges

**BECODISC B30C activated carbon stacked disc cartridges are ideal for the demanding liquid filtration applications in the chemical, pharmaceutical, cosmetics, food and beverage industries as well as in bioengineering.**

The specific advantages of BECODISC B30C activated carbon stacked disc cartridges:

- High adsorption power for decolorizing and aroma correcting.
- Complex fiber and cavity structure with a large interior surface for the widest range of applications and operating conditions.
- Easy to use and allow dust-free handling.

#### Adsorption through Activated Carbon

The activated carbon of BECODISC B30C is a micro-porous inert material, which is acid-washed and steam-activated. When products are cleaned or decolorized, a physical bond is created between the interior surfaces of the activated carbon and the unclean or colored substances. Since this bond is largely non-polar, there is a great affinity to organic molecules.

#### Factors Affecting the Adsorption Capacity

##### Filtration Speed

Adsorption processes are decisively affected by the contact time between the product and the adsorbing substance. The adsorption performance can thus be controlled by the speed of filtration. Slow filtration speeds and extended periods of contact result in optimum utilization of the adsorption capacity.



#### Water throughput BECO ACF 07



Conditions: Δ p = 14.5 psi (100 kPa, 1 bar), Medium: Water at 68 °F (20 °C)

#### Application Examples

- Decolorizing of chemicals and oils
- Decolorizing watery and alcoholic extracts
- Filtering of glucose, enzyme, vitamin, and antibiotic solutions
- Treating of cosmetics
- Taste and color correcting of beverages, spirits, and fruit juices

## Physical Data

This information is intended as a guideline for the selection of BECODISC stacked disc cartridges.

Type	Utilized BECO depth filter sheet	Thickness in (mm)	Ash content %	Bursting strength wet psi (kPa)	Water throughput at	
					$\Delta p = 14.5$ psi gpm/ft <sup>2</sup>	( $\Delta p = 100$ kPa* l/m <sup>2</sup> /min)
B30C	ACF 07	0.15 (3.8)	15.0	> 5.8 (40)	34.7	(1415)

The water throughput is a laboratory value characterizing the different BECO<sup>®</sup> depth filter sheets. It is not the recommended flow rate.

\* 100 kPa = 1 bar

## Chemical Data

Tested in accordance with § 177.2260 US Code of Federal Regulation of the US Food and Drug Administration (FDA).

Dry residue of the n-hexane extract 2-hour boiling time	< 0.5%
Dry residue of the extract with demineralized water 2-hour boiling time	< 4.0%
Dry residue of the extract with 5% acetic acid 2 hours at 194 °F (90 °C)	< 4.0%
Dry residue of the extract with 8% ethanol (v/v) 2 hours at 176 °F (80 °C)	< 4.0%
Dry residue of the extract with 50% ethanol (v/v) 2 hours at 140 °F (60 °C)	< 4.0%

## Components

BECODISC B30C activated stacked disc cartridges are manufactured from particularly pure materials, i.e. finely fibrillated cellulose fibers from deciduous and coniferous trees, cationic charge carriers, and high-quality diatomaceous earth.

## Recommendations for Avoiding Damage

BECODISC stacked disc cartridges can be used only in the specified flow direction. This applies to product filtering as well as sanitizing with hot water, and sterilizing with the stacked disc cartridges with saturated steam. In order to avoid damage to the filter cells, the system should be protected with a suitable non-return valve.

Refer to our insert included with each BECODISC stacked disc cartridge carton for detailed application information.

Depending on the filtered liquids, the operating temperature should not exceed 176 °F (80 °C). Please contact Eaton regarding filtration applications at higher temperatures.

## Intermediate Plates

If more than two BECODISC stacked disc cartridges (12" or 16") with double O-ring adapters are stacked in the housing, install a central spindle for safety reasons. In the event, more than one 16" BECODISC stacked disc cartridge (flat adapter/double O-ring adapter) is used in the housing, Eaton recommends the installation of stainless steel intermediate plates between the BECODISC stacked disc cartridges.

## Sanitizing and Sterilizing (Optional)

### Sterilizing with Hot Water

The hot water temperature should be 185 °F (85 °C) to 203 °F (95 °C). A differential pressure of 21.8 psi (150 kPa, 1.5 bar) must not be exceeded when sterilizing with hot water.

Sterilization time: At least 30 minutes once a minimum temperature of 176 °F (80 °C) is reached at all filter openings. In the interest of energy conservation, the water may be circulated provided the specified temperatures are maintained.

### Sterilizing with Steam

The wetted BECODISC stacked disc cartridges can be sterilized with saturated steam up to a maximum temperature of **250 °F (121 °C)** as follows:

Steam quality:	The steam must be free of foreign particles and impurities.
Temperature:	Max. <b>250 °F (121 °C)</b> <b>(saturated steam)</b>
Duration:	Approx. 20 minutes after steam exits from all filter valves.
Rinsing:	After sterilizing with 1.23 gal/ft <sup>2</sup> (50 l/m <sup>2</sup> ) at 1.25 times the flow rate.

## Filter Preparation and Filtration

Unless already completed after sterilization, Eaton recommends pre-rinsing the closed filter with 1.23 gal/ft<sup>2</sup> (50 l/m<sup>2</sup>) of water at 1.25 times the flow rate prior to the first filtration. Depending on the application, this usually equals a rinsing time of 10 – 20 minutes. Test the entire filter for leakage at maximum operating pressure.

High-proof alcohol solutions and products that do not allow pre-rinsing with water should be circulated for 10 to 20 minutes. Dispose of the rinsing solution after rinsing.

### **Differential Pressure**

Terminate the filtration process once the maximum permitted differential pressure of 43.5 psi (300 kPa, 3 bar) is reached. A higher differential pressure could damage the depth filter sheet material. For safety reasons, a differential pressure of 21.8 psi (150 kPa, 1.5 bar) should not be exceeded in applications for separating microorganisms.

### **Safety**

---

When used and handled correctly, there are no known unfavorable effects associated with this product.

Further safety information can be found in the relevant Material Safety Data Sheet, which can be downloaded from our website.

### **Waste Disposal**

---

Due to their composition, BECODISC stacked disc cartridges can be disposed of as harmless waste. Comply with relevant current regulations, depending on the filtered product.

### **Storage**

---

BECODISC stacked disc cartridges must be stored in a dry, odor-free, and well ventilated place.

Do not expose the BECODISC stacked disc cartridges to direct sunlight.

BECODISC stacked disc cartridges are intended for immediate use and should be used within 36 months after production date.

### **Available Formats**

---

BECODISC stacked disc cartridges are available with 12-inch and 16-inch diameters. Further information about filter areas and gasket types can be found in the current BECODISC stacked disc cartridge folder.

HS Customs Tariff: 84219900

### **Quality Assurance According to DIN EN ISO 9001**

---

Eaton's Begerow Product Line comprehensive Quality Management System has been certified according to DIN EN ISO 9001.

This certification verifies that a fully functioning comprehensive Quality Assurance System covering product development, contract controls, choice of suppliers, receiving inspections, production, final inspection, inventory management, and shipment has been implemented. Extensive quality assurance measures incorporate adherence to technical function criteria and chemical purity and quality recognized as safe under the German legislation governing the production of foods and beverages.

**North America**  
44 Apple Street  
Tinton Falls, NJ 07724  
Toll Free: 800 656-3344  
(North America only)  
Tel: +1 732 212-4700

**Europe/Africa/Middle East**  
Auf der Heide 2  
53947 Nettersheim, Germany  
Tel: +49 2486 809-0

Friedensstraße 41  
68804 Altlußheim, Germany  
Tel: +49 6205 2094-0

An den Nahewiesen 24  
55450 Langenlonsheim, Germany  
Tel: +49 6704 204-0

**China**  
No. 3, Lane 280,  
Linhong Road  
Changning District, 200335  
Shanghai, P.R. China  
Tel: +86 21 5200-0099

**Singapore**  
4 Loyang Lane #04-01/02  
Singapore 508914  
Tel: +65 6825-1668

**Brazil**  
Av. Julia Gaioli, 474 – Bonsucesso  
07251-500 – Guarulhos, Brazil  
Tel: +55 11 2465-8822

**For more information, please  
email us at [filtration@eaton.com](mailto:filtration@eaton.com)  
or visit [www.eaton.com/filtration](http://www.eaton.com/filtration)**

EN  
1 A 2.5.5.11  
08-2015

© 2015 Eaton. All rights reserved. All trademarks and registered trademarks are the property of their respective owners. All information and recommendations appearing in this brochure concerning the use of products described herein are based on tests believed to be reliable. However, it is the user's responsibility to determine the suitability for his own use of such products. Since the actual use by others is beyond our control, no guarantee, expressed or implied, is made by Eaton as to the effects of such use or the results to be obtained. Eaton assumes no liability arising out of the use by others of such products. Nor is the information herein to be construed as absolutely complete, since additional information may be necessary or desirable when particular or exceptional conditions or circumstances exist or because of applicable laws or government regulations.